

Your

Needle Guide


Intramuscular • Subcutaneous



Kills Phobia /Anxiety

Easy and Comfortable









Don't like needles / Don't like to inject

Self-Injection News... Needle fear and self-injection struggling Is eliminated......!

- Precision medicine dose delivered
- Outpatient programs delivere system
- No more anxiety and needle phobia
- Right injection angle (90 degree)
- Correct speed force
- Ensures stability at needle penetration
- LED light to prevent over-pressing
- Controlled tissue expansion
- Less injection pain

- Less side effects
- Controlled Delivery Target Deposit
- Aspiration check/avoid direct deposit in vessel
- Needle shield eliminates needle exposure
- Safety lock and Childproof
- Personalized IM/SC Injections
- Keeps needle steady
- Your dose is safe here
- Easy To Use

Injector for convenience and safety

Quality of care for Quality of life

Specifications & Features

90° Reusable NeedleGuide Intramuscular/Subcutaneous Dimensions 146 x 39 mm Syringes range 05 – 5ml Needle 30G and larger Needle range 8 – 40mm

